

System Requirements: Sensei Cloud

Sensei Cloud Workstation

- Desktop/Laptop Computer
 - Intel core i3 or later
 - 8GB RAM recommended
 - 1Gb of free space recommended
 - Monitor capable of 1440x900 screen resolution (1920x1080 recommended)
- Browsers:
 - **Windows:** Google Chrome or MS Edge web browser (latest version recommended)
 - **MacOS:** Google Chrome.
(note: Safari and Firefox can be used after selecting 'Continue' on the Browser Warning, but are not recommended currently)
- Internet Connection
 - 10 Mbps down / 10 Mbps up (25 Mbps down / 10 Mbps up recommended)
 - 4G mobile broadband (5G/LTE recommended)
- Supported Operating Systems
 - Actively supported version of Windows 10 or higher
- Imaging Software
 - CS Imaging Version 8 (Windows only) with Imaging Suite 3.3.0.4
See CSI System Requirements at <https://www.carestreamdental.com/en-us/resource-library/>.
- Scanners
 - Sensei Supports Twain compatible Scanners

Sensei Cloud Workstation with Cloud Imaging

- Desktop/Laptop Computer
 - Intel core i3 or later
 - 8GB RAM recommended
 - 8Gb of free space recommended
 - Monitor capable of 1440 x 900 screen resolution (1920x1080 recommended)
- Google Chrome web browser (latest version recommended)
- Internet Connection
 - 10 Mbps down / 10 Mbps up minimum
 - 4G mobile broadband (5G/LTE recommended)
- Supported Operating Systems
 - Actively supported version of Windows 10 or higher

Integration Requirements

If using an integrated product or service with Sensei Cloud, you must meet the system requirements of that product or service in addition to the Sensei Cloud System Requirements (this document).

Imaging Software

- CS Imaging Version 8 (Windows only) with Imaging Suite 3.3.0.4

See CSI System Requirements at <https://www.carestreamdental.com/en-us/resource-library/>.